

NANOGP8 Antibody



PACO01136

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

NANOGP8

Uniprot

Q6NSW7

Synonyms:

NANOGP8; Putative homeobox protein NANOGP8


Immunogen:

Synthesized peptide derived from the Internal region of human Nanog P8.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using Nanog P8 Polyclonal Antibody.