

NEGR1 Antibody



PACO01152

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

NEGR1

Uniprot

Q7Z3B1

Synonyms:

NEGR1; IGLON4; Neuronal growth regulator 1; IgLON family member 4


Immunogen:

Synthesized peptide derived from the Internal region of human NEGR1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using NEGR1 Polyclonal Antibody.