

# TSPAN9 Antibody



PACO01156

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

TSPAN9

**Uniprot**

O75954

**Synonyms:**

TSPAN9; NET5; Tetraspanin-9; Tspan-9; Tetraspan NET-5

**Immunogen:**


Synthesized peptide derived from the Internal region of human NET-5.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of A549 cells using NET-5 Polyclonal Antibody.