

NXPH3 Antibody



PACO01160

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

NXPH3

Uniprot

O95157

Synonyms:

NXPH3; KIAA1159; NPH3; Neurexophilin-3


Immunogen:

Synthesized peptide derived from the Internal region of human Neurexophilin-3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HT29 cells using Neurexophilin-3 Polyclonal Antibody.