

NKX2-4 Antibody



PACO01176

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

NKX2-4

Uniprot

Q9H2Z4

Synonyms:

NKX2-4; NKX2D; Homeobox protein Nkx-2.4; Homeobox protein NK-2 homolog D


Immunogen:

Synthesized peptide derived from the N-terminal region of human Nkx-2.4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using Nkx-2.4 Polyclonal Antibody.