

# OR10J1 Antibody



PACO01207

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

OR10J1

**Uniprot**

P30954/Q2M1M8/Q5VSV1

**Synonyms:**

Olfactory receptor 10J1

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 10J1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of Jurkat cells using Olfactory receptor 10J1 Polyclonal Antibody.

