

OR1A1 Antibody



PACO01213

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

OR1A1

Uniprot

Q9P1Q5

Synonyms:

OR1A1; Olfactory receptor 1A1; Olfactory receptor 17-7; OR17-7; Olfactory receptor OR17-11


Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 1A1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using Olfactory receptor 1A1 Polyclonal Antibody.