

OR1B1 Antibody



PACO01214

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR1B1

Uniprot

Q8NGR6

Synonyms:

OR1B1; Olfactory receptor 1B1; Olfactory receptor 9-B; OR9-B; Olfactory receptor OR9-26


Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 1B1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using Olfactory receptor 1B1 Polyclonal Antibody.