

# OR1D5/OR1D4 Antibody



PACO01215

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

OR1D5/OR1D4

**Uniprot**

P58170/P47884

**Synonyms:**

OR1D5; Olfactory receptor 1D5; Olfactory receptor 17-31; OR17-31; OR1D4; Olfactory receptor 1D4; Olfactory receptor 17-30; OR17-30

**Immunogen:**

Synthesized peptide derived from the Internal region of human Olfactory receptor 1D4/1D5.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of COLO205 cells using Olfactory receptor 1D4/1D5 Polyclonal Antibody.

