

OR1L6 Antibody



PACO01216

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR1L6

Uniprot

Q8NGR2

Synonyms:

OR1L6; OR1L7; Olfactory receptor 1L6; Olfactory receptor 1L7; Olfactory receptor OR9-30

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 1L6.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HeLa cells using Olfactory receptor 1L6 Polyclonal Antibody.

