

OR2AG1/OR2AG2 Antibody



PACO01219

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR2AG1/OR2AG2

Uniprot

Q9H205/A6NM03

Synonyms:

OR2AG1; OR2AG3; Olfactory receptor 2AG1; HT3; Olfactory receptor 2AG3; Olfactory receptor OR11-79; OR2AG2; OR2AG2P; Olfactory receptor 2AG2

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 2AG1/2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using Olfactory receptor 2AG1/2 Polyclonal Antibody.