

OR2T2/OR2T35 Antibody



PACO01226

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR2T2/OR2T35

Uniprot

Q6IF00/Q8NGX2

Synonyms:

OR2T2; OR2T2P; Olfactory receptor 2T2; Olfactory receptor OR1-43; OR2T35; Olfactory receptor 2T35; Olfactory receptor OR1-66

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 2T2/35.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HepG2 cells using Olfactory receptor 2T2/35 Polyclonal Antibody.

