

OR4P4 Antibody



PACO01236

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

OR4P4

Uniprot

Q8NGL7

Synonyms:

OR4P4; OR4P3P; Olfactory receptor 4P4; Olfactory receptor 4P3

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 4P4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HeLa cells using Olfactory receptor 4P4 Polyclonal Antibody.

