

OR9Q2 Antibody



PACO01264

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR9Q2

Uniprot

Q8NGE9

Synonyms:

OR9Q2; OR9Q2P; Olfactory receptor 9Q2

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 9Q2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HT29 cells using Olfactory receptor 9Q2 Polyclonal Antibody.