

# OMG Antibody



PACO01265

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

OMG

**Uniprot**

P23515

**Synonyms:**

OMG; OMGP; Oligodendrocyte-myelin glycoprotein

**Immunogen:**

Synthesized peptide derived from the Internal region of human OMG.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of 293 cells using OMG Polyclonal Antibody.

