

# TNFRSF11B Antibody



PACO01267

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300,  
IF:1:200-1:1000

**Protein Background:****Gene ID:**

TNFRSF11B

**Uniprot**

O00300

**Synonyms:**

TNFRSF11B; OCIF; OPG; Tumor necrosis factor receptor superfamily member 11B;  
Osteoclastogenesis inhibitory factor; Osteoprotegerin

**Immunogen:**

Synthesized peptide derived from the N-terminal region of human OPG.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of 293 cells using OPG Polyclonal Antibody.

