

PXN Antibody



PACO01305

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

PXN

Uniprot

P49023

Synonyms:

PXN; Paxillin


Immunogen:

Synthesized peptide derived from human Paxillin around the non-phosphorylation site of Y31.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COS7 cells using Paxillin Polyclonal Antibody.