

OLFML2B Antibody



PACO01332

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

OLFML2B

Uniprot

Q68BL8

Synonyms:

OLFML2B; Olfactomedin-like protein 2B; Photomedin-2


Immunogen:

Synthesized peptide derived from the Internal region of human Photomedin-2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc cells using Photomedin-2 Polyclonal Antibody.