

PIK3R1 Antibody



PACO01334

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

PIK3R1

Uniprot

P27986

Synonyms:

PIK3R1; GRB1; Phosphatidylinositol 3-kinase regulatory subunit alpha; PI3-kinase regulatory subunit alpha; PI3K regulatory subunit alpha; PtdIns-3-kinase regulatory subunit alpha; Phosphatidylinositol 3-kinase 85 kDa regulatory subunit alph


Immunogen:

Synthesized peptide derived from human PI 3-kinase p85alpha around the non-phosphorylation site of Y607.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 453 293 AD293 HELA HepG2 cells using PI 3-kinase p85alpha Polyclonal Antibody.