

# PIGX Antibody



PACO01340

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

PIGX

**Uniprot**

Q8TBF5

**Synonyms:**

PIGX; Phosphatidylinositol-glycan biosynthesis class X protein; PIG-X

**Immunogen:**


Synthesized peptide derived from the C-terminal region of human PIG-X.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of Jurkat cells using PIG-X Polyclonal Antibody.