

PPP4R1L Antibody



PACO01360

Product Information

Size:

50ug

Reactivity:

Human, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

PPP4R1L

Uniprot

Q9P1A2

Synonyms:

PPP4R1L; C20orf192; Putative serine/threonine-protein phosphatase 4 regulatory subunit 1-like

Immunogen:

Synthesized peptide derived from the Internal region of human PPP4R1L.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of COS7 cells using PPP4R1L Polyclonal Antibody.

