

AKT1S1 Antibody



PACO01362

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

AKT1S1

Uniprot

Q96B36

Synonyms:

AKT1S1; PRAS40; Proline-rich AKT1 substrate 1; 40 kDa proline-rich AKT substrate

Immunogen:

Synthesized peptide derived from human PRAS40 around the non-phosphorylation site of T246.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of 293 cells using PRAS40 Polyclonal Antibody.

