

PRSS33 Antibody



PACO01368

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

PRSS33

Uniprot

Q8NF86

Synonyms:

PRSS33; Serine protease 33; Serine protease EOS

Immunogen:

Synthesized peptide derived from the Internal region of human PRSS33.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of Jurkat cells using PRSS33 Polyclonal Antibody.

