

# BBC3 Antibody



PACO01380

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:****Gene ID:**

BBC3

**Uniprot**

Q9BXH1

**Synonyms:**

BBC3; PUMA; Bcl-2-binding component 3; JFY-1; p53 up-regulated modulator of apoptosis

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human PUMA.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of 293 cells using PUMA Polyclonal Antibody.

