

RAB29 Antibody



PACO01386

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

RAB29

Uniprot

O14966

Synonyms:

RAB7L1; Ras-related protein Rab-7L1; Rab-7-like protein 1

Immunogen:

Synthesized peptide derived from the Internal region of human Rab 7L1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HT29 cells using Rab 7L1 Polyclonal Antibody.

