

RACGAP1 Antibody



PACO01391

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

RACGAP1

Uniprot

Q9H0H5

Synonyms:

RACGAP1; KIAA1478; MGCRCACGAP; Rac GTPase-activating protein 1; Male germ cell RacGap; MgcRacGAP; Protein CYK4 homolog; CYK4; HsCYK-4

Immunogen:

Synthesized peptide derived from human Rac GAP1 around the non-phosphorylation site of S387.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of COS7 cells using Rac GAP1 Polyclonal Antibody.

