

# RRAD Antibody



PACO01393

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300,  
IF:1:200-1:1000

**Protein Background:****Gene ID:**

RRAD

**Uniprot**

P55042

**Synonyms:**

RRAD; RAD; GTP-binding protein RAD; RAD1; Ras associated with diabetes

**Immunogen:**


Synthesized peptide derived from the Internal region of human Rad GTPase.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of MCF7 cells using Rad GTPase Polyclonal Antibody.