

RAD50 Antibody



PACO01395

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

RAD50

Uniprot

Q92878

Synonyms:

RAD50; DNA repair protein RAD50; hRAD50

Immunogen:

Synthesized peptide derived from the Internal region of human Rad50.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HuvEc cells using Rad50 Polyclonal Antibody.

