

# RFPL4A Antibody



PACO01423

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

RFPL4A

**Uniprot**

A6NLU0

**Synonyms:**

RFPL4A; RFPL4; RNF210; Ret finger protein-like 4A; RING finger protein 210

**Immunogen:**


Synthesized peptide derived from the Internal region of human RFPL4A.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of HuvEc cells using RFPL4A Polyclonal Antibody.