

RFX3 Antibody



PACO01424

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

RFX3

Uniprot

P48380

Synonyms:

RFX3; Transcription factor RFX3; Regulatory factor X 3


Immunogen:

Synthesized peptide derived from the C-terminal region of human RFX3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of RAW264.7 cells using RFX3 Polyclonal Antibody.