

RPL10L Antibody



PACO01435

Product Information

Size:

50ug

Reactivity:

Human, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

RPL10L

Uniprot

Q96L21

Synonyms:

RPL10L; 60S ribosomal protein L10-like


Immunogen:

Synthesized peptide derived from the C-terminal region of human Ribosomal Protein L10L.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using Ribosomal Protein L10L Polyclonal Antibody.