

RPL39 Antibody



PACO01446

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

RPL39

Uniprot

P62891

Synonyms:

RPL39; 60S ribosomal protein L39


Immunogen:

Synthesized peptide derived from the N-terminal region of human Ribosomal Protein L39.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HepG2 cells using Ribosomal Protein L39 Polyclonal Antibody.