

SERPINB12 Antibody



PACO01490

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

SERPINB12

Uniprot

Q96P63

Synonyms:

SERPINB12; Serpin B12


Immunogen:

Synthesized peptide derived from the Internal region of human Serpin B12.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HeLa cells using Serpin B12 Polyclonal Antibody.