

SOX8/SOX9/SOX17/SOX18 Antibody



PACO01526

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

SOX8/SOX9/SOX17/SOX18

Uniprot

P57073/P48436/Q9H6I2/P35713

Synonyms:

SOX8; Transcription factor SOX-8; SOX9; Transcription factor SOX-9; SOX17; Transcription factor SOX-17; SOX18; Transcription factor SOX-18


Immunogen:

Synthesized peptide derived from the Internal region of human Sox-8/9/17/18.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COS7 cells using Sox-8/9/17/18 Polyclonal Antibody.