

SPIC Antibody



PACO01529

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

SPIC

Uniprot

Q8N5J4

Synonyms:

SPIC; Transcription factor Spi-C

Immunogen:

Synthesized peptide derived from the Internal region of human Spi-C.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HeLa cells using Spi-C Polyclonal Antibody.

