

SYT11 Antibody



PACO01564

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

SYT11

Uniprot

Q9BT88

Synonyms:

SYT11; KIAA0080; Synaptotagmin-11; Synaptotagmin XI; SytXI

Immunogen:

Synthesized peptide derived from the Internal region of human Synaptotagmin XI.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of CoLo205 cells using Synaptotagmin XI Polyclonal Antibody.