

SNCA Antibody



PACO01569

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

SNCA

Uniprot

P37840

Synonyms:

SNCA; NACP; PARK1; Alpha-synuclein; Non-A beta component of AD amyloid; Non-A4 component of amyloid precursor; NACP

Immunogen:

Synthesized peptide derived from human Synuclein-alpha around the non-phosphorylation site of S129.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HuvEc cells using Synuclein-alpha Polyclonal Antibody.

