

TAS2R39 Antibody



PACO01574

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

TAS2R39

Uniprot

P59534

Synonyms:

TAS2R39; Taste receptor type 2 member 39; T2R39; Taste receptor type 2 member 57; T2R57

Immunogen:

Synthesized peptide derived from the Internal region of human T2R39.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HepG2 cells using T2R39 Polyclonal Antibody.

