

# TAS2R8 Antibody



PACO01575

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

TAS2R8

**Uniprot**

Q9NYW2

**Synonyms:**

TAS2R8; Taste receptor type 2 member 8; T2R8; Taste receptor family B member 5; TRB5

**Immunogen:**

Synthesized peptide derived from the Internal region of human T2R8.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of COLO205 cells using T2R8 Polyclonal Antibody.

