

TACC3 Antibody



PACO01579

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

TACC3

Uniprot

Q9Y6A5

Synonyms:

TACC3; ERIC1; Transforming acidic coiled-coil-containing protein 3; ERIC-1

Immunogen:

Synthesized peptide derived from the C-terminal region of human TACC3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HuvEc cells using TACC3 Polyclonal Antibody.

