

NXF1 Antibody



PACO01585

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

NXF1

Uniprot

Q9UBU9

Synonyms:

NXF1; TAP; Nuclear RNA export factor 1; Tip-associated protein; Tip-associating protein; mRNA export factor TAP

Immunogen:

Synthesized peptide derived from the N-terminal region of human TAP.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HepG2 cells using TAP Polyclonal Antibody.

