

TBC1D3/TBC1D3B/TBC1D3C Antibody



PACO01588

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

TBC1D3/TBC1D3B/TBC1D3C

Uniprot

Q8IZP1/A6NDS4/Q6IPX1

Synonyms:

TBC1D3; PRC17; TBC1D3A; TBC1 domain family member 3; Prostate cancer gene 17 protein; Protein TRE17-alpha; Rab GTPase-activating protein PRC17; TBC1D3B; TBC1D3I; TBC1 domain family member 3B; TBC1D3C; TBC1D3D; TBC1 domain family member 3C

Immunogen:

Synthesized peptide derived from the C-terminal region of human TBC1D3A/B/C.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HeLa cells using TBC1D3A/B/C Polyclonal Antibody.

