

TBX1 Antibody



PACO01589

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

TBX1

Uniprot

O43435

Synonyms:

TBX1; T-box transcription factor TBX1; T-box protein 1; Testis-specific T-box protein

Immunogen:

Synthesized peptide derived from the C-terminal region of human TBX1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of COLO205 cells using TBX1 Polyclonal Antibody.

