

TBX10 Antibody



PACO01590

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

TBX10

Uniprot

O75333

Synonyms:

TBX10; TBX7; T-box transcription factor TBX10; T-box protein 10

Immunogen:

Synthesized peptide derived from the Internal region of human TBX10.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of HT29 cells using TBX10 Polyclonal Antibody.

