

# TBX3/TBX2 Antibody



PACO01593

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

TBX3/TBX2

**Uniprot**

O15119/Q13207

**Synonyms:**

TBX3; T-box transcription factor TBX3; T-box protein 3; TBX2; T-box transcription factor TBX2; T-box protein 2

**Immunogen:**

Synthesized peptide derived from the Internal region of human TBX2/3.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of A549 cells using TBX2/3 Polyclonal Antibody.

