

TCEAL1 Antibody



PACO01594

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IHC:1:200-1:1000

Protein Background:

Gene ID:

TCEAL1

Uniprot

Q15170

Synonyms:

TCEAL1; SIIR; Transcription elongation factor A protein-like 1; TCEA-like protein 1; Nuclear phosphoprotein p21/SIIR; Transcription elongation factor S-II protein-like 1

Immunogen:

Synthesized peptide derived from the Internal region of human TCEAL1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using TCEAL1 Polyclonal Antibody.