

# TEF Antibody



PACO01603

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

TEF

**Uniprot**

Q10587

**Synonyms:**

TEF; KIAA1655; Thyrotroph embryonic factor

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human TEF.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

Western Blot analysis of mouse brain cells using TEF Polyclonal Antibody.

