

THAP4 Antibody



PACO01622

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

THAP4

Uniprot

Q8WY91

Synonyms:

THAP4; CGI-36; PP238; THAP domain-containing protein 4

Immunogen:

Synthesized peptide derived from the Internal region of human THAP4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

Western Blot analysis of COS7 cells using THAP4 Polyclonal Antibody.

