

TOP3B Antibody



PACO01645

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

TOP3B

Uniprot

O95985

Synonyms:

TOP3B; TOP3B1; DNA topoisomerase 3-beta-1; DNA topoisomerase III beta-1


Immunogen:

Synthesized peptide derived from the N-terminal region of human Topo IIIBeta-1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc cells using Topo II β -1 Polyclonal Antibody.