

XRCC4 Antibody



PACO01711

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

XRCC4

Uniprot

Q13426

Synonyms:

XRCC4; DNA repair protein XRCC4; X-ray repair cross-complementing protein 4


Immunogen:

Synthesized peptide derived from the C-terminal region of human XRCC4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of COLO205 cells using XRCC4 Polyclonal Antibody.